n o v o x ®

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COMPOSITION:

Ingredients: oxygen-enriched olive oil (98%)
Excipients: natural pear fragrance (2%)
Substrate: polyurethane and polyester
The device may contain traces (<0.4%) of the following allergens: CITRAL, GERANIOL, CITRONELLOL, LINALOOL, LIMONENE

PRODUCT DESCRIPTION

NovoX® CUP mini is a medical device for the treatment of surgical, traumatic or ulcerative lesions and breast scars.

By creating a micro-environment which favours the activation of microcirculation but which is at the same time unfavourable to pathogen proliferation, and thanks also to the protective, barrier and soothing action of oxygen-enriched olive oil, NovoX® CUP aids the healing and correct cicatrization of lesions.

NovoX® CUP is supplied in the form of a single-use cup coated on the inside with an oily gel. The shape of the NovoX® CUP is specifically designed to fit the anatomical form of the breast.

NovoX® CUP is available in different sizes to suit different breast sizes.
For size selection, please refer to the indications on the package.

PACKAGE CONTENTS

2-piece box: 2 individually packaged impregnated substrates; 1 instruction manual (this leaflet)
6-piece box: 6 individually packaged impregnated substrates; 1 instruction manual (this leaflet)

INTENDED USES

NovoX® CUP is intended for the treatment of:

- Surgical wounds after oncological surgery (total or partial) and reconstructive breast surgery
- Cosmetic breast surgery wounds
- Hypertrophic scars
- Ulcers due to neoplastic diseases
- Post-radiotherapy ulcers and lesions
- Non-healing surgical wounds
- Superficial, medium and deep burns
- Traumatic wounds (e.g. abrasions, cuts and lacerations)
- Superficial lesions (including fissures)

NovoX® CUP is a non-sterile medical device. However, any previous contamination of the lesion is not a contraindication to the use of the device as tests have shown the intrinsic ability of the oxygen-enriched oil to create an environment unfavourable to the proliferation of pathogens commonly found on skin lesions, and is therefore inherently protected against contamination by such pathogens.

CONTRAINDICATIONS

NovoX® CUP must not be used on patients with known hypersensitivity or who have had an allergic reaction to the device or its components.

NovoX® CUP must not be used in case of previous serious skin reactions to wound dressings.

In case of known hypersensitivity to dressings, NovoX® CUP must not be used further to medical assessment.

NovoX® CUP is not recommended when complete coverage of the lesion is not possible.

PREGNANCY AND BREASTFEEDING

NovoX® CUP is a device for topical use not specifically tested for use during pregnancy. If you are or might be pregnant, we recommend that you do not use the device without medical advice.

NovoX® CUP is a device for topical use not specifically tested for use during breastfeeding. We recommend that

you do not use the device without medical advice.

INTENDED USERS

The device may be used by non-professional users.

In the case of extensive lesions or lesions involving the deep layers of the skin, the health care personnel should provide prior explanation on how to apply the dressing correctly.

Follow the instructions of the health care personnel regarding periodic check-ups.

SECONDARY DRESSINGS

To keep the device in position, the use of a bra is sufficient.

NovoX® CUP can be used with any bra as long as it is suitable for the needs of the post-operative period, following your doctor's instructions.

Follow the dressing protocol indicated by the health care personnel

DEVICE APPLICATION

Choose the most suitable NovoX® CUP size for the dimensions of your breast, following your doctor's advice. Make sure that the device covers the entire lesion area. For surgical wounds, apply NovoX® CUP immediately after the lesion is sutured.

Users should wash their hands thoroughly before and after use. Moreover, wear disposable gloves and use tools that are preferably sterile or otherwise subjected to careful disinfection when applying the product.

• Before applying the device, clean the lesion area as directed by health care personnel using physiological saline, sterile water or other appropriate solution. Pat dry with sterile gauze to remove excess water. Wipe off excess disinfectant if used.

• Open the package containing the device by pulling the flaps at the markings on the package, immediately before using the device.

• Take out the device, taking care not to touch the impregnated inside that will be in contact with the wounds.

• Fit the cup on the breast, taking care not to push the cup out of shape or crumple it excessively and check that the device covers the lesion or scar correctly.

• Put on your bra, taking care to keep the cup in position.
• Adjust the bra shoulder straps so that the cup adheres perfectly to the wounds or scars.

DEVICE REMOVAL

• To remove the device, simply lift it off the lesion. If removal is difficult, saturate with saline solution or sterile water and lift slowly.

DEVICE REPLACEMENT AND CONTINUED TREATMENT

The device should be replaced for the first time after 48 hours and then at least 2-3 times per week. The application of the device may be repeated for up to 30 days. However, the application frequency should be gradually reduced as the healing process is completed.

WARNINGS

Do not expose the device to intense heat or open flames. Do not use after the expiry date.

Do not reuse the device: the transmission of pathogens cannot be ruled out in the event of reuse.

Do not use the device if the package is opened or damaged.

Do not use the device if storage is not in accordance with these instructions

Do not use in conjunction with creams or medicinal products for topical use.

Do not use if other dressings have already been applied: this may reduce the effectiveness of the device and/or

cause local hypersensitivity reactions

Use on adhesive sutures may reduce the device's efficacy.

Use on anatomical locations other than the intended ones or on parts significantly different from the physiological anatomy, or excessive deformation or crumpling of the device, may reduce its efficacy and/or cause local hypersensitivity reactions.

A burning sensation or slight irritation immediately after application is normal and generally disappears within a short time. If the burning sensation is particularly intense or persistent, or if intense irritation appears, discontinue use and inform your doctor or paramedical staff.

The non-impregnated outside of the device is fairly water-repellent, but the device should not be wet or immersed in water.

Seepage of small amounts of oxygen-enriched oil from the edges of the dressing is normal.

Oxygen-enriched oil may stain fabrics

The pungent smell is an intrinsic characteristic, indicating the release of reactive oxygen species.

Keep out of the reach of children.

If the product is used on minors, adult supervision is recommended.

Promptly inform the manufacturer MOSS S.p.A. and the National Competent Authority if a serious accident occurs while using the device. Also inform the manufacturer MOSS S.p.A. of any side effects not mentioned in this leaflet.

If you require additional information, please contact the manufacturer, MOSS S.p.A.

DISPOSAL

Unused devices can be disposed of as unsorted waste.

The device and any secondary dressings after use should be disposed of as potentially infected biological materials.

STORAGE

Keep at a temperature below 25°C, in a cool, dry area far from sources of heat.

In case of high room temperatures the product should be kept in the refrigerator (2-8°C)

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 0477

Date of last revision of instruction leaflet:
Rev: MDR-0922 of 13/09/2022

 25°C	Keep at temperatures below 25 °C	REF	Product number	LOT	Lot number		Expiry date
MD	Medical device	UDI	Unique Device Identifier		Read the instructions for use prior to use		Do not reuse
	Read warnings prior to use		Do not use if the package is damaged		Keep away from light		Keep dry

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