NOVOX®

NOVOX®

NOVOX®

NOVOX®

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NOVOX®

Medical device for the treatment of lesions in correspondence of percutaneous access

EN

COMPOSITION:

Ingredients: oxygen-enriched olive oil (98%)
Excipients: natural pear fragrance (2%)
Substrate: polyurethane and polyester
The device may contain traces (<0.4 %) of the following allergens: CITRAL, GERANIOL, CITRONELLLOL, LINALOOL, LIMONENE.

PRODUCT DESCRIPTION

Novox® Spot is a medical device for the treatment of wounds, including ulcers, in correspondence of percutaneous accesses (e.g. tracheostomy, PEG, drainages, orthopedic fixators and similar, EXCLUDING VENOUS AND/OR ARTERIAL ACCESSSES AND INTESTINAL AND URINARY STOMA FOR COLLECTION BAGS).
By creating a microenvironment which favours activation of the microcirculation but is unfavourable for pathogen proliferation, and thanks also to the protective, barrier and soothing action of oxygen-enriched olive oil, NovoX® Spot aids the healing and correct cicatrization of wounds.

Novox® Spot is a medical device in the form of a disposable planar substrate impregnated on one side with an oily gel, with cutting and circular opening to adapt the device to the diameter of the percutaneous access.

PACKAGE CONTENTS

1-piece box: 1 individually packaged impregnated substrate; 1 instruction manual (this leaflet)
5-piece box: 5 individually packed impregnated substrates; 1 instruction manual (this leaflet)
10-piece box: 10 individually packaged impregnated substrates; 1 instruction manual (this leaflet)

DESTINATIONS OF USE

Novox® Spot is intended for the treatment of:
• Ulcers, or non-healing wounds lesions in correspondence of percutaneous accesses
• Surgical lesions in correspondence of percutaneous accesses

NovoX® Spot is a non-sterile medical device. However, any previous contamination of the lesion is not a contraindication to the use of the device, since tests have revealed the intrinsic capacity of the oxygen-enriched oil to create an environment unfavourable to the proliferation of the pathogens commonly found on skin lesions, making it inherently protected against contamination by these pathogens.

CONTRAINDICATIONS

Novox® Spot must not be used on patients with known hypersensitivity or who have had an allergic reaction to the device or its components.

Novox® Spot must not be used in case of previous serious skin reactions to the application of dressings. In case of known hypersensitivity to dressings, NovoX® Spot must be used further to medical assessment.
NovoX® Spot should not be used when complete coverage of the lesion is not possible.

PREGNANCY AND BREASTFEEDING

Novox® Spot is a topical device not specifically tested for use during pregnancy. If you are or might be pregnant, we recommend that you do not use the device without medical advice.

Novox® Spot is a topical device not specifically tested for use during breastfeeding. We recommend that you do not use the device without medical advice.

INTENDED USERS

The device can be used by non-professional users. In case of extensive lesions or lesions involving deep layers of the skin, it is recommended that healthcare personnel provide a prior explanation on how to apply the dressing correctly.
Follow the instructions of the healthcare personnel regarding the periodic check-ups.

SECONDARY DRESSINGS

To keep the device in place, secondary dressings may be used.

For the fixing of drainages or cannulas, secondary dressings can be used, taking care to ensure that the lumen of the tubes is not deformed or narrowed. It is recommended to follow in any case the dressing protocol indicated by the healthcare personnel.

DEVICE APPLICATION

Make sure that the device covers the entire lesion area. Users should wash their hands thoroughly before and after use. Moreover, wear disposable gloves and use tools that are preferably sterile or otherwise subjected to careful disinfection when applying the product.

• Before applying the device, clean the area of the lesion as instructed by the healthcare professional using physiological solution, sterile water or other appropriate solution. Dab with sterile gauze to remove excess water. If disinfectant is used, dry the excess disinfectant.

• Open the package containing the device by pulling the flaps at the stamped indication, immediately before using the device itself.

• Take out the device taking care not to touch the impregnated part that will come into contact with the wounds. Remove the protective film.

• Place the device on the site of the lesion, taking care not to deform or bend excessively and make sure that the device covers the entire lesion area.

• If necessary, it is possible to cut the device to reduce its size, after careful disinfection of the tools used (or using sterile tools) and taking care not to damage the support or the internal and external lining fabrics.

• Fix the device in place, where necessary using the secondary dressings indicated by the healthcare personnel and fix the drainage or cannula as instructed by the healthcare personnel.

REMOVING THE DEVICE

• To remove the device, it is sufficient to lift it from the injury. If removal is difficult, wet with saline solution or sterile water and lift slowly.

DEVICE REPLACEMENT AND CONTINUED TREATMENT

It is recommended to make the first replacement of the device after 48 hours and subsequent replacements at least 2-3 times a week. The application of the device can be repeated for a maximum time of 30 days. However, it is advisable to reduce the frequency of applications as you complete the healing process.

WARNINGS

Do not expose the device to sources of intense heat or open flames.

Do not use after the expiry date.

Do not reuse the device: in the case of reuse, transmission of pathogens cannot be excluded.

Do not use the device if the package is open or damaged. Do not use the device if storage is not in accordance with these instructions.

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0477
Date of last revision of product leaflet:
Rev: MDR-0922 of 13/09/2022

	Keep at temperatures below 25°C		Product number		Lot number		Expiry date - yyyy/mm/dd
	Medical device		Unique Device Identifier		Read the instructions for use prior to use		Do not reuse
	Read warnings prior to use		Do not use if the package is damaged		Keep away from light		Keep dry

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Para más información, se ruega contactar con el fabricante MOSS S.p.A.